



# F10

Fetal Doppler




## — A pocket-size fetal doppler that measures fetal heart rate and allows audio output of the fetal heart sound

- Each unit comes with a built-in speaker for direct sound output and an audio jack for earphones
- Each package includes a robust carrying case for easy transport and product protection
- Works with 2 standard AA batteries that allow long lasting battery operation (6 hours)
- Measures and indicates Fetal Heart Rate as early as 10 weeks into the pregnancy
- Highly sensitive Doppler probe, using a 2 MHz ultrasound frequency
- Clear and easy to read LCD display
- PC Interface Sound Card

## — F10 Ultrasound Doppler System

- Monitoring of fetal heart rate by health-care professionals and for personal use
- Physicians
- Pregnant Women
- Midwives
- Obstetricians





## TECHNICAL SPECIFICATION

Display	· LCD Back Light	· Low Battery Indication
	· Signal Strength Indication	
Heart Rate Range	· 50~240 bpm	
Accuracy	· ±2%	
Ultrasound Frequency	· 2MHz	
Ultrasound Intensity	· <10mW/cm <sup>2</sup>	
Speaker Output Power	· 1.2W (Built-in Loudspeaker)	
Audio Output for Earphone		
Auto Shut-off - Power	· 5 minutes	
Auto Shut-off - Sound	· 1 minute	
Power	· 1.5V Battery X 2 (Type: AA)	
Sensitivity	· 10~12 Weeks Onward	
PC Interface	· Sound Card	
Battery Life	· 360 minutes	
Accessories	· Robust Carrying Case	· Battery(1.5V AA × 2)
	· Gel(60ml)	

\* Specifications can be changed without prior notification.

- High-Sensitivity Doppler probe
- Ergonomic design
- Large LCD Display & Clear Sound
- User-friendly concept
- Continuous Wave Doppler
- Pulse Wave Doppler



Donghwa Medical Instrument Complex, 1650-1, Donghwa-ri,  
Munmak-eup, Wonju-si, Gangwon-do, 220-801 Korea

Tel : +82-33-742-5410 / Fax : +82-33-742-5418

Tel : +82-2-542-3375 / Fax : +82-2-542-7447

e-mail : info@mediana.co.kr (Overseas)

d.sales@mediana.co.kr (Domestic)

www.mediana.co.kr

